



Data sheet : **SONO MN E5**

Product description : Medical face mask with elastic bands and nickelfree nose clip on top of the filter.

Use : Medical face mask for the use by medical (or comparable) personal in contact with the patient. The mask reduces the probability of an infection of the patient and the user.
Not for the use during invasive operations (no splash resistance).

Physical properties

Outer layer (filter) : 100 % PP Meltblown, 100 g/qm
Inner layer (face) : 100 % PP Spunbond, 20 g/qm

Colour : white
Thickness : app. 1.0 mm
Bacterial filter efficiency : ≥ 98 % BFE (**Type II**)

Safety information

The mask is not sterile. All materials used are odourless and free of harmful substances.

Do not reuse! A reuse may cause severe or chronic infections.


Change before you visit a new patient.

Please change the mask if it is obviously damaged or moist.

Storage

Conditions : Please store in a dry and ventilated place at ambient temperature. Keep away from fire sources.

Max. storage time : 36 month



Issued 02 November 2020

You will find this material data sheet on our homepage: www.dfabielefeld.de & www.dfa-gmbh.com