



**Data sheet** : **SONO MN E6**

**Product description** : Medical face mask with elastic bands and nickelfree nose clip on top of the filter.

**Use** : Medical face mask for the use by medical (or comparable) personal in contact with the patient. The mask reduces the probability of an infection of the patient and the user including splash resistance.

### Physical properties

Outer layer (anti-splash) : 100 % PP Spunbond, 40 g/qm  
Outer layer (filter) : 100 % PP Meltblown, 80 g/qm  
Inner layer (face) : 100 % PP Spunbond, 20 g/qm

Colour : Blue and white  
Thickness : app. 1.3 mm  
Bacterial filter efficiency :  $\geq 98$  % BFE (**Type IIR**)

### Safety information

The mask is not sterile. All materials used are odourless and free of harmful substances.

Do not reuse! A reuse may cause severe or chronic infections.

Change before you visit a new patient.

Please change the mask if it is obviously damaged or moist.

### Storage

Conditions : Please store in a dry and ventilated place at ambient temperature. Keep away from fire sources.

Max. storage time : 36 month

Issued October 12, 2020

You will find this material data sheet on our homepage: [www.dfabielefeld.de](http://www.dfabielefeld.de) & [www.dfa-gmbh.com](http://www.dfa-gmbh.com)